wearable breast pump ZOMEE FIT



Now you can have it all – a Zomee **hospital-grade strength** motor in a wearable pump! Thoughtfully designed for quick and effective expression with Zomee's hands-free wearable pumps. Rechargeable for hours of use and easy-to-use technology that doesn't require an app.

Inside The Box

- Two (2) 25mm Breast Shields
- Two (2) 21mm Breast Shield Inserts
- Two (2) 17mm Breast Shield Inserts
- Two (2) Pump Motors
- Two (2) Breast Shield Bodies
- Two (2) Diaphragms
- Two (2) 150ml Collection Containers
- Two (2) Collection Container Lids
- Four (4) Seal Rings
- Four (4) Duckbill Valves
- One (2) Type-C Charging Cable
- Two (2) Bra Adjusters
- One (1) Collection Container Brush
- One (1) English + Spanish Manual
- One (1) Nipple Ruler + Guide

Features

- Hospital Grade Strength
- Leak Proof
- 2+ Hours Battery Life (Quick Charging Protocol -Using PD/QC charger to charge 10mins, the pump could work for 30mins)
- Reaches 260MMHG- 26~ 100CPM
- 9 Levels Stimulation
- 9 Levels Expression
- 9 Levels 2-Phase Mode (mimicks way baby feeds from mom giving mom extra stimulation for quicker let down)
- Memory Function (Remembers the last setting)
- Auto LED Screen Lock (After 30 seconds, screen goes off automatically)
- 2 Mins Auto-Switch from Stimulation to Expression
- Display LED Screen (display of mode, battery status & levels)
- Noise level- ≤48dB
- Soft Liquid Silicone Breast Shield- sizes 17/21/25MM in the box
- Anti-Backflow System
- BPA + Phthalate Free
- Easy to Clean and Assemble



Worn inside your bra. Silent. Hands-Free. Comfortable.